

## PRINTED RED TABARDS for Nursing & Care Home Staff

Machine Washable Red Tabards printed with “Do Not Disturb. Nurse On Medication Round”.

Helping to identify nurses on medication rounds.  
These red tabards are manufactured in the UK from a performance fabric making them an excellent choice for NHS & Care home nursing staff.

Made from a quality PVC/elastomer, bonded to a knitted polyester backing.  
It has an easily wipeable surface and can also be Machine Washed at 40°C.  
Bound with black binding with 25mm wide cotton tape side ties.

Printed on the front & back with a choice of:

“Do Not Disturb. Nurse On Medication Round” or “Do Not Disturb. Nurse On Drug Round”.

### Washing Instructions








Wash separately with Standard bio washing powder or liquid upto 40°C.

Surface can be wiped with an antibacterial spray.

Do not Tumble-dry

Do not Dry clean

Do not use solvent-based products on the printed surface as this may result in the print ink coming off.

|  |                    |   |
|--|--------------------|---|
|    | Machine Wash, Warm | Initial water temperature should not exceed 40C or 105F.  |
|  | Do Not Dryclean    | Garment may not be commercially drycleaned.   |
|  | Do Not Iron        | Item may not be smoothed or finished with an iron.  |
|  | Do Not Bleach      | No bleach product may be used. The garment is not colorfast or structurally able to withstand any bleach. |
|  | Do Not Wring       | Do not wring.   |
|  | Do Not Tumble Dry  | A machine dryer may not be used.  |
|  | Drip Dry           | Hang dripping wet garment from line or bar, in or out doors, without hand shaping or smoothing            |



Unit 313 Jubilee Trade Centre  
130 Pershore Street  
Birmingham  
B5 6ND

Tel: 0121 666 6422

Email: [sales@kovamanufacturing.co.uk](mailto:sales@kovamanufacturing.co.uk)

Web: [www.kovamanufacturing.co.uk](http://www.kovamanufacturing.co.uk)

Web: [www.logopens.co.uk/kova](http://www.logopens.co.uk/kova)